

TENT5C 抗原（重组蛋白）

中文名称：TENT5C 抗原（重组蛋白）

英文名称：TENT5C Antigen (Recombinant Protein)

别名：FAM46C

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 192-391 amino acids of human TENT5C

技术规格：

Full name:	terminal nucleotidyltransferase 5C
Synonyms:	FAM46C
Swissprot:	Q5VWP2
Gene Accession:	BC036516
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Probable nucleotidyltransferase that may act as a non-canonical poly(A) RNA polymerase. (Microbial infection) Seems to enhance replication of some viruses, including yellow fever virus, in response to type I interferon.