

TEX11 抗原（重组蛋白）

中文名称：TEX11 抗原（重组蛋白）

英文名称：TEX11 Antigen (Recombinant Protein)

别名：testis expressed 11; TGC1; TSGA3; SPGFX2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 741-940 amino acids of human TEX11

技术规格：

Full name:	testis expressed 11
Synonyms:	TGC1; TSGA3; SPGFX2
Swissprot:	Q8IYF3
Gene Accession:	BC036016
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is X-linked and is expressed in only male germ cells. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.