

## TMEM9 抗原（重组蛋白）

中文名称：TMEM9 抗原（重组蛋白）

英文名称： TMEM9 Antigen (Recombinant Protein)

别名： DERM4; TMEM9A

储存： 冷冻（-20℃）

相关类别： 抗原

概述：

Fusion protein corresponding to a region derived from 21-183 amino acids of human TMEM9

技术规格：

<b>Full name:</b>	transmembrane protein 9
<b>Synonyms:</b>	DERM4; TMEM9A
<b>Swissprot:</b>	Q9P0T7
<b>Gene Accession:</b>	BC001106
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	May be involved in intracellular transport.