

TOR1AIP1 抗原（重组蛋白）

中文名称：TOR1AIP1 抗原（重组蛋白）

英文名称：TOR1AIP1 Antigen (Recombinant Protein)

别名：LAP1; LAP1B; LGMD2Y

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 384-583 amino acids of human TOR1AIP1

技术规格：

Full name:	torsin 1A interacting protein 1
Synonyms:	LAP1; LAP1B; LGMD2Y
Swissprot:	Q5JTV8
Gene Accession:	BC023247
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a type 2 integral membrane protein that binds A- and B-type lamins. The encoded protein localizes to the inner nuclear membrane and may be involved in maintaining the attachment of the nuclear membrane to the nuclear lamina during cell division. Alternate splicing results in multiple transcript variants.