

TOR1A 抗原（重组蛋白）

中文名称：TOR1A 抗原（重组蛋白）

英文名称：TOR1A Antigen (Recombinant Protein)

别名：torsin family 1 member A; DQ2; DYT1

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 21-197 amino acids of human TOR1A

技术规格：

Full name:	torsin family 1 member A
Synonyms:	DQ2; DYT1
Swissprot:	O14656
Gene Accession:	BC014484
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the AAA family of adenosine triphosphatases (ATPases), is related to the Clp protease/heat shock family and is expressed prominently in the substantia nigra pars compacta. Mutations in this gene result in the autosomal dominant disorder, torsion dystonia 1.