

TRAPPC1 抗原（重组蛋白）

中文名称：TRAPPC1 抗原（重组蛋白）

英文名称：TRAPPC1 Antigen (Recombinant Protein)

别名：BET5; MUM2

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-145 amino acids of human TRAPPC1

技术规格：

Full name:	trafficking protein particle complex 1
Synonyms:	BET5; MUM2
Swissprot:	Q9Y5R8
Gene Accession:	BC032717
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene product plays a role in vesicular transport of proteins to the Golgi apparatus from the endoplasmic reticulum. The encoded protein is a component of the multisubunit transport protein particle (TRAPP) complex. Alternative splicing results in multiple transcript variants.