

## TRO 抗原（重组蛋白）

中文名称：TRO 抗原（重组蛋白）

英文名称：TRO Antigen (Recombinant Protein)

别名：MAGED3; MAGE-d3

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-200 amino acids of human TRO

技术规格：

<b>Full name:</b>	trophinin
<b>Synonyms:</b>	MAGED3; MAGE-d3
<b>Swissprot:</b>	Q12816
<b>Gene Accession:</b>	BC026914
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a membrane protein that mediates cell adhesion between trophoblastic cells and the epithelial cells of the endometrium. The encoded protein participates in cell signaling during embryo implantation, and may also be involved in cancer formation. This gene is located near several other closely related genes on chromosome X. Alternative splicing results in multiple transcript variants.