

兔抗 MAGEB10 多克隆抗体

中文名称：兔抗 MAGEB10 多克隆抗体

英文名称： Anti-MAGEB10 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

抗 原： MAGEB10

宿 主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the B subfamily of the melanoma associated antigen protein family. The encoded protein is specifically expressed in testis and tumor cells.
Applications:	ELISA, IHC
Name of antibody:	MAGEB10
Immunogen:	Fusion protein of human MAGEB10
Full name:	melanoma antigen family B, 10
SwissProt:	Q96LZ2
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human thyroid cancer and human brain
IHC Recommend dilution:	100-300

