

兔抗 GPC6 多克隆抗体

中文名称：兔抗 GPC6 多克隆抗体

英文名称：Anti-GPC6 rabbit polyclonal antibody

别名：OMIMD1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GPC6

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The glypicans comprise a family of glycosylphosphatidylinositol-anchored heparan sulfate proteoglycans, and they have been implicated in the control of cell growth and cell division. The glypican encoded by this gene is a putative cell surface coreceptor for growth factors, extracellular matrix proteins, proteases and anti-proteases.
Applications:	ELISA, IHC
Name of antibody:	GPC6
Immunogen:	Fusion protein of human GPC6
Full name:	glypican 6

Synonyms :	OMIMD1
SwissProt:	Q9Y625
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain
IHC Recommend dilution:	50-100

