

## 兔抗 GPKOW 多克隆抗体

中文名称：兔抗 GPKOW 多克隆抗体

英文名称： Anti-GPKOW rabbit polyclonal antibody

别名： G-patch domain and KOW motifs; T54; Spp2; GPATC5; GPATCH5

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： GPKOW

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a putative RNA-binding protein containing G-patch and KOW (Kyprides, Ouzounis, Woese) domains. The encoded protein interacts directly with protein kinase A and protein kinase X and is also found associated with the spliceosome.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	GPKOW
<b>Immunogen:</b>	Fusion protein of human GPKOW
<b>Full name:</b>	G-patch domain and KOW motifs
<b>Synonyms:</b>	T54; Spp2; GPATC5; GPATCH5

<b>SwissProt:</b>	Q92917
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	100-500
<b>WB Predicted band size:</b>	53 kDa
<b>WB Positive control:</b>	293T and Jurkat cell lysates
<b>WB Recommended dilution:</b>	1000-5000



