

兔抗 GPR1 多克隆抗体

中文名称: 兔抗 GPR1 多克隆抗体

英文名称: Anti-GPR1 rabbit polyclonal antibody

别名: G protein-coupled receptor 1

抗原: GPR1

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Receptor for the inflammation-associated leukocyte chemoattractant chemerin/RARRES2 suggesting a role for this receptor in the regulation of inflammation. Can act as a coreceptor for HIV-1.
Applications:	ELISA, WB, IHC
Name of antibody:	GPR1
Immunogen:	Synthetic peptide of human GPR1
Full name:	G protein-coupled receptor 1
SwissProt:	P46091
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and Human thyroid cancer

IHC Recommend dilution:	20-100
WB Predicted band size:	41 kDa
WB Positive control:	Hela cell lysate
WB Recommended dilution:	200-1000



