

## 兔抗 GPR15 多克隆抗体

中文名称：兔抗 GPR15 多克隆抗体

英文名称：Anti-GPR15 rabbit polyclonal antibody

别名：BOB

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GPR15

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a G protein-coupled receptor that acts as a chemokine receptor for human immunodeficiency virus type 1 and 2. The encoded protein localizes to the cell membrane.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	GPR15
<b>Immunogen:</b>	Synthetic peptide of humanGPR15
<b>Full name:</b>	G protein-coupled receptor 15
<b>Synonyms :</b>	BOB
<b>SwissProt:</b>	P49685

<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human esophagus cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	10-50
<b>WB Predicted band size:</b>	41 kDa
<b>WB Positive control:</b>	HepG2, 231, Hela and Lovo cells
<b>WB Recommended dilution:</b>	500-2000



