

## 兔抗 GPR17 多克隆抗体

中文名称: 兔抗 GPR17 多克隆抗体

英文名称: Anti-GPR17 rabbit polyclonal antibody

别名: G protein-coupled receptor 17

储存: 冷冻 (-20℃)

抗原: GPR17

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Dual specificity receptor for uracil nucleotides and cysteinyl leukotrienes (CysLTs). Signals through G(i) and inhibition of adenylyl cyclase. May mediate brain damage by nucleotides and CysLTs following ischemia.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GPR17
<b>Immunogen:</b>	Synthetic peptide of human GPR17
<b>Full name:</b>	G protein-coupled receptor 17
<b>SwissProt:</b>	Q13304
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer

IHC Recommend dilution:

30-150



