

## 兔抗 GPR149 多克隆抗体

中文名称: 兔抗 GPR149 多克隆抗体

英文名称: Anti-GPR149 rabbit polyclonal antibody

别名: G protein-coupled receptor 149; IEDA; PGR10

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: GPR149

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Orphan receptor.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	GPR149
<b>Immunogen:</b>	Synthetic peptide of human GPR149
<b>Full name:</b>	G protein-coupled receptor 149
<b>Synonyms:</b>	IEDA; PGR10
<b>SwissProt:</b>	Q86SP6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	81 kDa
<b>WB Positive control:</b>	Human prostate tissue lysate

WB Recommended dilution:

200-1000



