

ZNF688 抗原（重组蛋白）

中文名称：ZNF688 抗原（重组蛋白）

英文名称：ZNF688 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述：

Fusion protein corresponding to a region derived from 1-276 amino acids of human ZNF688

技术规格：

Full name:	zinc finger protein 688
Swissprot:	P0C7X2
Gene Accession:	BC018997
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May be involved in transcriptional regulation.