

## 兔抗 GPRC6A 多克隆抗体

中文名称: 兔抗 GPRC6A 多克隆抗体

英文名称: Anti-GPRC6A rabbit polyclonal antibody

别名: G protein-coupled receptor class C group 6 member A; GPCR; bA86F4.3

抗原: GPRC6A

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Members of family C of the G protein-coupled receptor (GPCR) superfamily, such as GPRC6A, are characterized by an evolutionarily conserved amino acid-sensing motif linked to an intramembranous 7-transmembrane loop region. Several members of GPCR family C, including GPRC6A, also have a long N-terminal domain (summary by Pi et al., 2005 [PubMed 16199532]).
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	GPRC6A
<b>Immunogen:</b>	Synthetic peptide of human GPRC6A
<b>Full name:</b>	G protein-coupled receptor class C group 6 member A

<b>Synonyms:</b>	GPCR; bA86F4.3
<b>SwissProt:</b>	Q5T6X5
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	20-100
<b>WB Predicted band size:</b>	105 kDa
<b>WB Positive control:</b>	Hela cell lysate
<b>WB Recommended dilution:</b>	200-1000



