

兔抗 GPRIN1 多克隆抗体

中文名称：兔抗 GPRIN1 多克隆抗体

英文名称： Anti-GPRIN1 rabbit polyclonal antibody

别名： GRIN1

抗原： GPRIN1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	May be involved in neurite outgrowth By similarity.
Applications:	WB
Name of antibody:	GPRIN1
Immunogen:	Synthesized peptide derived from internal of human GPRIN1.
Full name:	G protein regulated inducer of neurite outgrowth 1
Synonyms :	GRIN1
SwissProt:	Q7Z2K8
WB Predicted band size:	102 kDa
WB Positive control:	HT-29 cells lysate
WB Recommended dilution:	500-3000

