

兔抗 GREB1 多克隆抗体

中文名称: 兔抗 GREB1 多克隆抗体

英文名称: Anti-GREB1 rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃)

宿 主: Rabbit

抗 原: GREB1

反应种属: Human, Mouse

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene is an estrogen-responsive gene that is an early response gene in the estrogen receptor-regulated pathway. It is thought to play an important role in hormone-responsive tissues and cancer. Three alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	GREB1
Immunogen:	Synthetic peptide of human GREB1
Full name:	growth regulation by estrogen in breast cancer 1
SwissProt:	Q4ZG55
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human brain

IHC Recommend dilution:	50-200
-------------------------	--------

