

兔抗 GSC 多克隆抗体

中文名称: 兔抗 GSC 多克隆抗体

英文名称: Anti-GSC rabbit polyclonal antibody

别名: gooseoid homeobox; SAMS

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: GSC

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the bicoid subfamily of the paired (PRD) homeobox family of proteins. The encoded protein acts as a transcription factor and may be autoregulatory. A similar protein in mice plays a role in craniofacial and rib cage development during embryogenesis.
Applications:	ELISA, IHC
Name of antibody:	GSC
Immunogen:	Synthetic peptide of human GSC
Full name:	gooseoid homeobox
Synonyms:	SAMS
SwissProt:	P56915

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and Human esophagus cancer
IHC Recommend dilution:	20-100



