

兔抗 GSK3B 多克隆抗体

中文名称：兔抗 GSK3B 多克隆抗体

英文名称：Anti-GSK3B rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：GSK3B

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	The protein encoded by this gene is a serine-threonine kinase, belonging to the glycogen synthase kinase subfamily. It is involved in energy metabolism, neuronal cell development, and body pattern formation. Polymorphisms in this gene have been implicated in modifying risk of Parkinson disease, and studies in mice show that overexpression of this gene may be relevant to the pathogenesis of Alzheimer disease. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	GSK3B
Immunogen:	Synthetic peptide of human GSK3B

Full name:	glycogen synthase kinase 3 beta
SwissProt:	P49841
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100

