

## 兔抗 GPR63 多克隆抗体

- 中文名称：兔抗 GPR63 多克隆抗体
- 英文名称：Anti-GPR63 rabbit polyclonal antibody
- 别名：G protein-coupled receptor 63; PSP24B; PSP24(beta)
- 相关类别：一抗
- 储存：冷冻 (-20℃)
- 宿主：Rabbit
- 抗原：GPR63
- 反应种属：Human, Mouse
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a G protein-coupled receptor. Multiple alternatively spliced variants, encoding the same protein, have been identified. Orphan receptor. May play a role in brain function.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GPR63
<b>Immunogen:</b>	Synthetic peptide of human GPR63
<b>Full name:</b>	G protein-coupled receptor 63
<b>Synonyms:</b>	PSP24B; PSP24(beta)
<b>SwissProt:</b>	Q9BZJ6
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human cervical cancer and Human colorectal cancer
<b>IHC Recommend dilution:</b>	25-100



