

## 兔抗 GPR75 多克隆抗体

中文名称：兔抗 GPR75 多克隆抗体

英文名称：Anti-GPR75 rabbit polyclonal antibody

别名：G protein-coupled receptor 75; GPRchr2; WI31133

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：GPR75

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	GPR75 is a member of the G protein-coupled receptor family. GPRs are cell surface receptors that activate guanine-nucleotide binding proteins upon the binding of a ligand.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	GPR75
<b>Immunogen:</b>	Synthetic peptide of human GPR75
<b>Full name:</b>	G protein-coupled receptor 75
<b>Synonyms:</b>	GPRchr2; WI31133
<b>SwissProt:</b>	O95800
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	50-300

