

兔抗 MAP3K12 多克隆抗体

中文名称：兔抗 MAP3K12 多克隆抗体

英文名称： Anti-MAP3K12 rabbit polyclonal antibody

别名： DLK; MUK; ZPK; ZPKP1; MEKK12

抗原： MAP3K12

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a member of the serine/threonine protein kinase family. This kinase contains a leucine-zipper domain and is predominately expressed in neuronal cells. The phosphorylation state of this kinase in synaptic terminals was shown to be regulated by membrane depolarization via calcineurin. This kinase forms heterodimers with leucine zipper containing transcription factors, such as cAMP responsive element binding protein (CREB) and MYC, and thus may play a regulatory role in PKA or retinoic acid induced neuronal differentiation. Alternatively spliced transcript variants encoding different proteins have been described.

Applications:	ELISA, IHC
Name of antibody:	MAP3K12
Immunogen:	Fusion protein of human MAP3K12
Full name:	mitogen-activated protein kinase kinase kinase 12
Synonyms :	DLK; MUK; ZPK; ZPKP1; MEKK12
SwissProt:	Q12852
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and human thyroid cancer
IHC Recommend dilution:	25-100

