

兔抗 MAP3K7(Ab-184) 多克隆抗体

中文名称：兔抗 MAP3K7(Ab-184) 多克隆抗体

英文名称：Anti-MAP3K7(Ab-184) rabbit polyclonal antibody

别名：TAK1; MEKK7; TGF1a

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：MAP3K7(Ab-184)

应种属：Human Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the serine/threonine protein kinase family. This kinase mediates the signaling transduction induced by TGF beta and morphogenetic protein (BMP), and controls a variety of cell functions including transcription regulation and apoptosis. In response to IL-1, this protein forms a kinase complex including TRAF6, MAP3K7P1/TAB1 and MAP3K7P2/TAB2; this complex is required for the activation of nuclear factor kappa B.
Applications:	WB
Name of antibody:	MAP3K7(Ab-184)
Immunogen:	Synthesized non-phosphopeptide derived from human TAK

	1 around the phosphorylation site of threonine 184 (I-Q-T(p)-H-M).
Full name:	mitogen-activated protein kinase kinase kinase 7
Synonyms :	TAK1; MEKK7; TGF1a
SwissProt:	O43318
WB Predicted band size:	67 kDa
WB Positive control:	HepG2 cells lysate
WB Recommended dilution:	500-3000

