

兔抗 HCRT/pan Orexin 多克隆抗体

中文名称：兔抗 HCRT/pan Orexin 多克隆抗体

英文名称：Anti-HCRT rabbit polyclonal antibody

别名：OX; PPOX; NRCLP1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HCRT/pan Orexin

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a hypothalamic neuropeptide precursor protein that gives rise to two mature neuropeptides, orexin A and orexin B, by proteolytic processing. Orexin A and orexin B, which bind to orphan G-protein coupled receptors HCRTR1 and HCRT R2, function in the regulation of sleep and arousal. This neuropeptide arrangement may also play a role in feeding behavior, metabolism, and homeostasis.
Applications:	ELISA, IHC
Name of antibody:	HCRT/pan Orexin
Immunogen:	Synthetic peptide of human HCRT/pan Orexin
Full name:	hypocretin (orexin) neuropeptide precursor

Synonyms:	OX; PPOX; NRCLP1
SwissProt:	O43612
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	25-100

