

ABCB8 抗原（重组蛋白）

中文名称： ABCB8 抗原（重组蛋白）

英文名称： ABCB8 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 394-693 amino acids of human ABCB8

技术规格

Full name:	ATP-binding cassette sub-family B member 8
Synonyms:	MABC1, M-ABC1, EST328128
Swissprot:	Q9NUT2
Gene Accession:	BC141836
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance as well as antigen presentation. The function of this half-transporter has not yet been determined; however, it may involve the compartmentalization and transport of heme, as well as peptides, from the mitochondria to the nucleus.

s and cytosol. This protein may also play a role in the transport of phospholipids into mitochondrial membranes.