

兔抗 HCRTR2 多克隆抗体

中文名称：兔抗 HCRTR2 多克隆抗体

英文名称：Anti-HCRTR2 rabbit polyclonal antibody

别名：OX2R

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HCRTR2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a G-protein coupled receptor involved in the regulation of feeding behavior. The encoded protein binds the hypothalamic neuropeptides orexin A and orexin B. A related gene (HCRTR1) encodes a G-protein coupled receptor that selectively binds orexin A.
Applications:	ELISA, IHC
Name of antibody:	HCRTR2
Immunogen:	Synthetic peptide of human HCRTR2
Full name:	hypocretin (orexin) receptor 2
Synonyms :	OX2R
SwissProt:	O43614

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human lung cancer and Human ovarian cancer
IHC Recommend dilution:	50-200

