

兔抗 MAP4K1 多克隆抗体

中文名称：兔抗 MAP4K1 多克隆抗体

英文名称： Anti-MAP4K1 rabbit polyclonal antibody

别名： HPK1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MAP4K1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene belongs to the protein kinase superfamily, STE Ser/Thr protein kinase family, STE20 subfamily. It plays a role in the response to environmental stress and hematopoietic lineage decisions and growth regulation. Appears to act upstream of the JUN N-terminal pathway.
Applications:	ELISA, IHC
Name of antibody:	MAP4K1
Immunogen:	Synthetic peptide of human MAP4K1
Full name:	mitogen-activated protein kinase kinase kinase 1
Synonyms :	HPK1
SwissProt:	Q92918

ELISA Recommended dilution:	1000-5000
IHC positive control:	Human brain and Human gastric cancer
IHC Recommend dilution:	25-100

