

兔抗 MAP4K2 多克隆抗体

中文名称：兔抗 MAP4K2 多克隆抗体

英文名称： Anti-MAP4K2 rabbit polyclonal antibody

别名： GCK; BL44; RAB8IP

相关类别： 一抗

抗原： MAP4K2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the serine/threonine protein kinase family. Although this kinase is found in many tissues, its expression in lymphoid follicles is restricted to the cells of germinal centre, where it may participate in B-cell differentiation. This kinase can be activated by TNF-alpha, and has been shown to specifically activate MAP kinases. This kinase is also found to interact with TNF receptor-associated factor 2 (TRAF2), which is involved in the activation of MAP3K1/MEKK1.
Applications:	ELISA, WB
Name of antibody:	MAP4K2

Immunogen:	Fusion protein of human MAP4K2
Full name:	mitogen-activated protein kinase kinase kinase kinase 2
Synonyms :	GCK; BL44; RAB8IP
SwissProt:	Q12851
ELISA Recommended dilution:	2000-10000
WB Predicted band size:	92 kDa
WB Positive control:	Human fetal brain tissue
WB Recommended dilution:	1000-5000

