

兔抗 MAP4K4 多克隆抗体

中文名称：兔抗 MAP4K4 多克隆抗体

英文名称： Anti-MAP4K4 rabbit polyclonal antibody

别名： HGK; NIK; MEKKK4; FLH21957; HEL-S-31

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： MAP4K4

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene is a member of the serine/threonine protein kinase family. This kinase has been shown to specifically activate MAPK8/JNK. The activation of MAPK8 by this kinase is found to be inhibited by the dominant-negative mutants of MAP3K7/TAK1, MAP2K4/MKK4, and MAP2K7/MKK7, which suggests that this kinase may function through the MAP3K7-MAP2K4-MAP2K7 kinase cascade, and mediate the TNF-alpha signaling pathway.
Applications:	WB, IF
Name of antibody:	MAP4K4
Immunogen:	Synthesized peptide derived from internal of human MEKKK

	4.
Full name:	mitogen-activated protein kinase kinase kinase kinase 4
Synonyms :	HGK; NIK; MEKKK4; FLH21957; HEL-S-31
SwissProt:	O95819
WB Predicted band size:	142 kDa
WB Positive control:	COLO cells and Huvec cells lysates
WB Recommended dilution:	500-3000
IF positive control:	A549 cells
IF Recommend dilution:	100-500



