

## 兔抗 MAP4K5 多克隆抗体

中文名称：兔抗 MAP4K5 多克隆抗体

英文名称： Anti-MAP4K5 rabbit polyclonal antibody

别名： KHS; GCKR; KHS1; MAPKKK5

抗原： MAP4K5

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the serine/threonine protein kinase family, that is highly similar to yeast SPS1/STE20 kinase. Yeast SPS1/STE20 functions near the beginning of the MAP kinase signal cascades that is essential for yeast pheromone response. This kinase was shown to activate Jun kinase in mammalian cells, which suggested a role in stress response. Two alternatively spliced transcript variants encoding the same protein have been described for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MAP4K5

<b>Immunogen:</b>	Fusion protein of human MAP4K5
<b>Full name:</b>	mitogen-activated protein kinase kinase kinase kinase
<b>Synonyms :</b>	KHS; GCKR; KHS1; MAPKKKK5
<b>SwissProt:</b>	Q9Y4K4
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human ovarian cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	15-50

