

ANXA5 抗原（重组蛋白）

中文名称：ANXA5 抗原（重组蛋白）

英文名称： ANXA5 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 49-263 amino acids of human ANXA5

技术规格

Full name:	annexin A5
Synonyms:	PP4, ANX5, ENX2, RPRGL
Swissprot:	P08758
Gene Accession:	BC018671
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene belongs to the annexin family of calcium-dependent phospholipid binding proteins some of which have been implicated in membrane-related events along exocytotic and endocytotic pathways. Annexin 5 is a phospholipase A2 and protein kinase C inhibitory protein with calcium channel activity and a potential role in cellular signal transduction, inflammation, growth and differentiation. Annexin 5 has also been described as placental anticoagulant protein I, vascular anticoagulant-alpha, endonexin II, lipocortin V, placental protein 4 and anchorin CII. The gene spans 29 kb containing 13 exons, and encodes a single transcript of approximate

ly 1.6 kb and a protein product with a molecular weight of about 3
5 kDa.