

兔抗 MAPK1 多克隆抗体

中文名称：兔抗 MAPK1 多克隆抗体

英文名称： Anti-MAPK1 rabbit polyclonal antibody

别名： ERK; p38; p40; p41; ERK2; ERT1; ERK-2; MAPK2; PRKM1; PRKM2; P42MAPK; p41mapk; p42-MAPK

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： MAPK1

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the MAP kinase family. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. The activation of this kinase requires its phosphorylation by upstream kinases.
Applications:	WB
Name of antibody:	MAPK1
Immunogen:	Synthesized peptide derived from internal of human p42

	MAPK.
Full name:	mitogen-activated protein kinase 1
Synonyms :	ERK; p38; p40; p41; ERK2; ERT1; ERK-2; MAPK2; PRKM1; PRKM2; P42MAPK; p41mapk; p42-MAPK
SwissProt:	P28482
WB Predicted band size:	41 kDa
WB Positive control:	Jurkat cells and 293 cells lysates
WB Recommended dilution:	500-3000

