

兔抗 HDAC10 多克隆抗体

中文名称：兔抗 HDAC10 多克隆抗体

英文名称：Anti-HDAC10 rabbit polyclonal antibody

抗原：HDAC10

别名：HD10

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the histone deacetylase family, members of which deacetylate lysine residues on the N-terminal part of the core histones. Histone deacetylation modulates chromatin structure, and plays an important role in transcriptional regulation, cell cycle progression, and developmental events. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	HDAC10
Immunogen:	Fusion protein of human HDAC10
Full name:	histone deacetylase 10

Synonyms:	HD10
SwissProt:	Q969S8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	20-100



