

## 兔抗 HDC 多克隆抗体

中文名称：兔抗 HDC 多克隆抗体

英文名称：Anti-HDC rabbit polyclonal antibody

抗 原：HDC

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the group II decarboxylase family and forms a homodimer that converts L-histidine to histamine in a pyridoxal phosphate dependent manner. Histamine regulates several physiologic processes, including neurotransmission, gastric acid secretion, inflammation, and smooth muscle tone.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	HDC
<b>Immunogen:</b>	Fusion protein of human HDC
<b>Full name:</b>	histidine decarboxylase
<b>SwissProt:</b>	P19113
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human colon cancer and Human ovarian cancer

IHC Recommend dilution:	5-20
-------------------------	------

