

兔抗 MAPK7 多克隆抗体

中文名称：兔抗 MAPK7 多克隆抗体

英文名称：Anti-MAPK7 rabbit polyclonal antibody

别名：BMK1; ERK4; ERK5; PRKM7

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MAPK7

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the MAP kinase family. MAP kinases act as an integration point for multiple biochemical signals, and are involved in a wide variety of cellular processes such as proliferation, differentiation, transcription regulation and development. This kinase is specifically activated by mitogen-activated protein kinase kinase 5 (MAP2K5/MEK5). It is involved in the downstream signaling processes of various receptor molecules including receptor tyrosine kinases, and G protein-coupled receptors.
Applications:	ELISA, WB, IHC
Name of antibody:	MAPK7

Immunogen:	Fusion protein of human MAPK7
Full name:	mitogen-activated protein kinase 7
Synonyms :	BMK1; ERK4; ERK5; PRKM7
SwissProt:	Q13164
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and human lung cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	88 kDa
WB Positive control:	Human fetal muscle tissue
WB Recommended dilution:	200-1000



