

## 兔抗 HDHD2 多克隆抗体

中文名称：兔抗 HDHD2 多克隆抗体

英文名称：Anti-HDHD2 rabbit polyclonal antibody

别名：HEL-S-301; 3110052N05Rik

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HDHD2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	HDHD2 (haloacid dehalogenase-like hydrolase domain containing 2) is also known as DKFZp564D1378 and is a 259 amino acid protein that is expressed as two isoforms produced by alternative splicing. HDHD2 belongs to the HAD-like hydrolase superfamily, which contains a group of hydrolase enzymes
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	HDHD2
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	haloacid dehalogenase-like hydrolase domain containing 2
<b>Synonyms :</b>	HEL-S-301; 3110052N05Rik
<b>SwissProt:</b>	Q9H0R4

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and human cervical cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	29 kDa
<b>WB Positive control:</b>	HEPG2 cell lysate
<b>WB Recommended dilution:</b>	200-1000





