

## 兔抗 HDHD5 多克隆抗体

中文名称: 兔抗 HDHD5 多克隆抗体

英文名称: Anti-HDHD5 rabbit polyclonal antibody

别名: haloacid dehalogenase like hydrolase domain containi

抗原: HDHD5

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Candidate gene for the Cat Eye Syndrome (CES), a developmental disorder associated with the duplication of a 2 Mb region of 22q 11.2. Duplication usually takes in the form of a supernumerary bi satellited isodicentric chromosome, resulting in four copies of the region (represents an inv dup(22)(q11)). CES is characterized clinically by the combination of coloboma of the iris and anal atresia with fistula, downslanting palpebral fissures, preauricular tags and/or pits, frequent occurrence of heart and renal malformations, and normal or near-normal mental development.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	HDHD5

<b>Immunogen:</b>	Synthetic peptide of human HDHD5
<b>Full name:</b>	haloacid dehalogenase like hydrolase domain containing 5
<b>Synonyms:</b>	CECR5
<b>SwissProt:</b>	Q9BXW7
<b>WB Predicted band size:</b>	46 kDa
<b>WB Positive control:</b>	Jurkat,HEPG2 and Hela cell,Human testis tissue,231 cell lysates
<b>WB Recommended dilution:</b>	500-2000

