

兔抗 HEATR1 多克隆抗体

中文名称: 兔抗 HEATR1 多克隆抗体

英文名称: Anti-HEATR1 rabbit polyclonal antibody

别名: HEAT repeat containing 1; BAP28; UTP10

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: HEATR1

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Ribosome biogenesis factor. Involved in nucleolar processing of pre-18S ribosomal RNA. Required for optimal pre-ribosomal RNA transcription by RNA polymerase I.
Applications:	ELISA, IHC
Name of antibody:	HEATR1
Immunogen:	Fusion protein of human HEATR1
Full name:	HEAT repeat containing 1
Synonyms:	BAP28; UTP10
SwissProt:	Q9H583
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human prostate cancer
IHC Recommend dilution:	50-100

