

## 兔抗 HEATR9 多克隆抗体

中文名称：兔抗 HEATR9 多克隆抗体

英文名称：Anti-HEATR9 rabbit polyclonal antibody

别名：HEAT repeat containing 9; C17orf66

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HEATR9

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The function of this protein remains unknown.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	HEATR9
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	HEAT repeat containing 9
<b>Synonyms:</b>	C17orf66
<b>SwissProt:</b>	A2RTY3
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	50-300

