

BCL2L13 抗原（重组蛋白）

中文名称： BCL2L13 抗原（重组蛋白）

英文名称： BCL2L13 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 57-356 amino acids of human BCL2L13

技术规格

Full name:	BCL2-like 13 (apoptosis facilitator)
Synonyms:	MIL1; BCL-RAMBO; Bcl2-L-13
Swissprot:	Q9BXK5
Gene Accession:	BC003032
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a mitochondrially-localized protein with conserved B-cell lymphoma 2 homology motifs. Overexpression of the encoded protein results in apoptosis. Alternatively spliced transcript variants have been observed for this gene.