

兔抗 MAPRE3 多克隆抗体

中文名称：兔抗 MAPRE3 多克隆抗体

英文名称： Anti-MAPRE3 rabbit polyclonal antibody

别名： EB3; RP3; EBF3; EBF3-S

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： MAPRE3

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the RP/EB family of genes. The protein localizes to the cytoplasmic microtubule network and binds APCL, a homolog of the adenomatous polyposis coli tumor suppressor gene.
Applications:	ELISA, WB, IHC
Name of antibody:	MAPRE3
Immunogen:	Fusion protein of human MAPRE3
Full name:	microtubule-associated protein, RP/EB family, member 3
Synonyms :	EB3; RP3; EBF3; EBF3-S
SwissProt:	Q9UPY8
ELISA Recommended dilution	2000-5000

:	
IHC positive control:	Human breast cancer and human prostate cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	32 kDa
WB Positive control:	Human fetal brain and mouse brain tissue
WB Recommended dilution:	500-2000



