

兔抗 HEBP2 多克隆抗体

中文名称:	兔抗 HEBP2 多克隆抗体
英文名称:	Anti-HEBP2 rabbit polyclonal antibody
别名:	heme binding protein 2; PP23; SOUL; C6orf34; C6ORF34B
相关类别:	一抗
储存:	冷冻 (-20℃)
宿主:	Rabbit
抗原:	HEBP2
反应种属:	Human, Mouse
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is found predominately in the cytoplasm, where it plays a role in the collapse of mitochondrial membrane potential (MMP) prior to necrotic cell death. The encoded protein enhances outer and inner mitochondrial membrane permeabilization, especially under conditions of oxidative stress.
Applications:	ELISA, IHC
Name of antibody:	HEBP2
Immunogen:	Fusion protein of human HEBP2
Full name:	heme binding protein 2
Synonyms:	PP23; SOUL; C6orf34; C6ORF34B
SwissProt:	Q9Y5Z4

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300

