

BIRC5 抗原（重组蛋白）

中文名称：BIRC5 抗原（重组蛋白）

英文名称： BIRC5 Antigen (Recombinant Protein)

别名： API4, EPR-1

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	baculoviral IAP repeat containing 5
Synonyms:	API4, EPR-1
Swissprot:	O15392
Gene Accession:	BC008718
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is a member of the inhibitor of apoptosis (IAP) gene family , which encode negative regulatory proteins that prevent apoptotic cell death. IAP family members usually contain multiple baculovirus IAP repeat (BIR) domains, but this gene encodes proteins with only a single BIR domain. The encoded proteins also lack a C-terminus RING finger domain. Gene expression is high during fetal development and in most tumors, yet low in adult tissues. Alternatively spliced transcript v

ariants encoding distinct isoforms have been found for this gene.