

兔抗 HECTD1 多克隆抗体

中文名称：兔抗 HECTD1 多克隆抗体

英文名称：Anti-HECTD1 rabbit polyclonal antibody

别名：HECT domain E3 ubiquitin protein ligase 1; EULIR

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HECTD1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Probable E3 ubiquitin-protein ligase which accepts ubiquitin from a n E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. May be required for development of the head mesenchyme and neural tube closure.
Applications:	ELISA, IHC
Name of antibody:	HECTD1
Immunogen:	Synthetic peptide of human HECTD1
Full name:	HECT domain E3 ubiquitin protein ligase 1
Synonyms:	EULIR
SwissProt:	Q9ULT8

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	10-50

