

兔抗 HECTD2 多克隆抗体

中文名称：兔抗 HECTD2 多克隆抗体

英文名称：Anti-HECTD2 rabbit polyclonal antibody

别名：HECT domain E3 ubiquitin protein ligase 2

相关类别：一抗

抗原：HECTD2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Probable E3 ubiquitin-protein ligase which accepts ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates.
Applications:	ELISA, IHC
Name of antibody:	HECTD2
Immunogen:	Synthetic peptide of human HECTD2
Full name:	HECT domain E3 ubiquitin protein ligase 2
SwissProt:	Q5U5R9
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human brain
IHC Recommend dilution:	50-100

