

## BTG2 抗原（重组蛋白）

中文名称：BTG2 抗原（重组蛋白）

英文名称： BTG2 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-158 amino acids of human BTG2

技术规格

<b>Full name:</b>	BTG anti-proliferation factor 2
<b>Synonyms:</b>	PC3; APRO1; TIS21
<b>Swissprot:</b>	P78543
<b>Gene Accession:</b>	BC105949
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a member of the BTG/Tob family. This family has structurally related proteins that appear to have antiproliferative properties. This encoded protein is involved in the regulation of the G1/S transition of the cell cycle.