

兔抗 HELLS 多克隆抗体

中文名称：兔抗 HELLS 多克隆抗体

英文名称：Anti-HELLS rabbit polyclonal antibody

别名：LSH; PASG; SMARCA6; Nbla10143

抗原：HELLS

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a lymphoid-specific helicase. Other helicases function in processes involving DNA strand separation, including replication, repair, recombination, and transcription. This protein is thought to be involved with cellular proliferation and may play a role in leukemogenesis. Alternatively spliced transcript variants encoding different isoforms have been identified.
Applications:	ELISA, WB, IHC
Name of antibody:	HELLS
Immunogen:	Fusion protein of human HELLS
Full name:	helicase, lymphoid-specific
Synonyms:	LSH; PASG; SMARCA6; Nbla10143

SwissProt:	Q9NRZ9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human esophagus cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	97 kDa
WB Positive control:	HEPG2 cell lysate
WB Recommended dilution:	500-2000



