

BTG4 抗原（重组蛋白）

中文名称：BTG4 抗原（重组蛋白）

英文名称： BTG4 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 27-206 amino acids of human BTG4

技术规格

Full name:	BTG family member 4
Synonyms:	PC3B
Swissprot:	Q9NY30
Gene Accession:	BC031045
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the BTG/Tob family. This family has structurally related proteins that appear to have antiproliferative properties. This encoded protein can induce G1 arrest in the cell cycle.