

## 兔抗 HELZ2 多克隆抗体

中文名称: 兔抗 HELZ2 多克隆抗体

英文名称: Anti-HELZ2 rabbit polyclonal antibody

别名: PDIP-1; PRIC285

抗原: HELZ2

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a nuclear transcriptional co-activator for peroxisome proliferator activated receptor alpha. The encoded protein contains a zinc finger and is a helicase that appears to be part of the peroxisome proliferator activated receptor alpha interacting complex. This gene is a member of the DNA2/NAM7 helicase gene family. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	HELZ2
<b>Immunogen:</b>	Synthesized peptide derived from internal of human PRIC285.
<b>Full name:</b>	helicase with zinc finger 2, transcriptional coactivator

<b>Synonyms :</b>	PDIP-1; PRIC285
<b>SwissProt:</b>	Q9BYK8
<b>WB Predicted band size:</b>	295 kDa
<b>WB Positive control:</b>	Jurkat cells lysate
<b>WB Recommended dilution:</b>	500-3000

